



Recombinant *Gossypium hirsutum* Apocytochrome f (petA)

| | |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-CF643522GHB |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | Q2L921 |
| Product Type | Transmembrane Protein |
| Immunogen Species | <i>Gossypium hirsutum</i> (Upland cotton) (<i>Gossypium mexicanum</i>) |
| Sequence | YPIFAQQGYENPREATGRIVCANCHLANKPVDIEVPQAVLPDTVFEAVVRIPYD MQLKQV LANGKKGALNVGAVLILPEGFELAPPDRISPENKKEKIGNLSFQNYRPTKKNILVI GPVPG KKYSEITFPILSPDPASNKDAHFLKYPIYVGGNRGRGQIYPDGNKSNNTVYNAT ATGIVS KIIRKEKGGYEITITDALDGHQVVDIIPPGPELLVSEGESIKLDQPLTINPNVGGF GQGD AEIVLQDPLRVQGLLFFLASIVFAQIFLVLKKKQFEKVQVSEMNF |
| Source | in vitro E.coli expression system |
| Target Names | petA |
| Protein Names | Recommended name: Apocytochrome f |
| Expression Region | 36-320 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |