



Recombinant *Saccharomyces cerevisiae* Golgi to ER traffic protein 2 (GET2)

Product Code	CSB-CF334663SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P40056
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MSELTEAEKRLLRERRQKKFSNGGASSRLNKITGQASSHLNAESPLDAPSAA KTTTPPAS VHSATPDIKEDSNVAPQLDLLKQLAAMQGQGTGKSTPQDSSTPDLLSLLSSMN TGMPSAE GTPSFGQAAPAAPINQAALDYHDYLLNRLKAWTILVKWVFFLLPYLYLITRPN SVWPAY AFTQSAWFAPLRNPSNFTRIFATFEFLSISIIYQLLKNVEHKSKIKNLQDTNKL V KLVSL VPEGVIPVANLKGKLITLLQYWDLLSMLITDISFVLIVLGLLTYL
Source	in vitro <i>E.coli</i> expression system
Target Names	GET2
Protein Names	Recommended name: Golgi to ER traffic protein 2 Alternative name(s): Guided entry of tail-anchored proteins 2 Hydroxyurea resistance protein 2 Required for meiotic nuclear division protein 7
Expression Region	1-285
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.