



Recombinant Escherichia coli O6:K15:H31 Uncharacterized protein yjiK (yjiK)

Product Code	CSB-CF606949EGY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q0T8X4
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Sequence	MTKSISLSKRIFVIVILFVIVAVCTFFVQSCARKSNHAASFQNYHATIDGKEIAGIT NNI SSLTWSAQSNLTFSTINKPAAIVEMTTNGDLIRTIPLDFVKDLETIEYIGDNQFVIS DER DYAIYVISLTPNSEVKILKKIKIPLQESPTNCGFEGLAYSRQDHTFWFFKEKNPIE VYKV NGLLSSNELHISKDKALQRQFTLDDVSGAEFNQQKNTLLVLSHESRALQEVTL VGEVIGE MSLTKGSRGLSHNIKQAEGVAMDASGNLYIVSEPNRFYRFTPQSSH
Source	in vitro E.coli expression system
Target Names	yjiK
Protein Names	Recommended name: Uncharacterized protein yjiK
Expression Region	1-286
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.