



Recombinant Uncharacterized protein yjiK (yjiK)

Product Code	CSB-CF813841EGX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8FA95
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Sequence	MTKSISLSKRIFVIVILFVIVAVCTFFVQSCARKSNHAASFQNYHATIDGKEIAGIT NNI SSLTWSAQSNLTFSTINKPAAIVEMTTNGDLIRTIPLDFVKDLETIEYIGDNQFVIS DER DYAIYVISLTPNSEVKILKKIKIPLQESPTNCGFEGLAYSRQDHTFWFFKEKNPIE VYKV NGLLSSNELHISKDKALQRQFTLDDVSGAEFNQQKNTLLVLSHESRALQEVTL VGEVIGE MSLTKGSRGLSHNIKQAEGVAMDASGNIYIVSEPNRFYRFTPQSSH
Source	in vitro E.coli expression system
Target Names	yjiK
Protein Names	Recommended name: Uncharacterized protein yjiK
Expression Region	1-286
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.