



Recombinant Arabidopsis thaliana RING-H2 finger protein ATL30 (ATL30)

Product Code	CSB-CF863767DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9FIR0
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MPIAKPINQNTTVPYPPQHYSKPPLVIILTVILLVVFFIGFFAIYFCKCFYHTLTEA WNH HYHNGLPENQIQAQQEPVQPPVNPGLPHIIQSYPPLFPFSSVKDLREDKYGLE CAICLLE FEEEHILLRLLTTCYHVFHQECIDQWLESNKTCPVCRRLDPNAPENIKELIIEVI QENA HENRDQEQTSTSNEVMLSRQSSGNNERKIETLPDKFSRSKTTGHSIVRNKPEE EDRYTLR LPDHVKIKVTRRHNNNQTESCISFGELVRNREGRFGEVSGQSLVPESGS
Source	in vitro E.coli expression system
Target Names	ATL30
Protein Names	Recommended name: RING-H2 finger protein ATL30
Expression Region	1-289
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.