



Recombinant Danio rerio Endoplasmic reticulum-Golgi intermediate compartment protein 1 (ergic1)

Product Code	CSB-CF703296DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q4V8Y6
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MSFDVRRFDIYRKVPKDLTQPTYTGAFISICCCVFMLFLFLSELTGFIATEIVNEL YVDD PDKDSGGKIDVSLNISLPNLHCDLVGLDIQDEMGRHEVGHIEISMKVPLNNGH GCRFEGE FSINKVPGNFHVSTHSATAQPQSPDMTHIIHKLAFGAKLQVQHVGAFNALGG ADRLQSN ALASHDYILKIVPTVYEELGGKQRFSYQYTVANKEYVAYSHTGRIIPAIWFRYDL SPITV KYTERRRPFYRFITTICAIIGGTFTVAGIIDSCIFTASEAWKKIQIGKMS
Source	in vitro E.coli expression system
Target Names	ergic1
Protein Names	Recommended name: Endoplasmic reticulum-Golgi intermediate compartment protein 1 Alternative name(s): ER-Golgi intermediate compartment 32 kDa protein Short name= ERGIC-32
Expression Region	1-290
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.