



Recombinant *Desulfotalea psychrophila* UPF0316 protein DP2912 (DP2912)

Product Code	CSB-CF731761DAAA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6AJ39
Product Type	Transmembrane Protein
Immunogen Species	<i>Desulfotalea psychrophila</i> (strain LSv54 / DSM 12343)
Sequence	MIMDFVSFFEPASLLYAAMVFLSRVVDVSFGTLRTISIVHGRTTMAFWLGFFEA AIWLVV VSTIVQTVTEQPLLIVYALGFATGNIVGIKFERMIALGHLVLRVISRHNSNEIAYE IRE QGYAVTAFKGEGRSGIVTELYIVCRRRDMKRILKKVMALDPKAFYVSERAGEV SNVYRPI MQPMTGWRRVLKKKYHSPSSALQGFSGEGSSPSSLSQLGRAERHFPGTVLR CRPYVPFLV ISFASEMRPPDTARAFSFGQRCGLARLTPPCISQRALHGTSLSGEALVSG
Source	in vitro E.coli expression system
Target Names	DP2912
Protein Names	Recommended name: UPF0316 protein DP2912
Expression Region	1-290
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.