



Recombinant Human Transmembrane O-methyltransferase (LRTOMT)

Product Code	CSB-CF837882HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8WZ04
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MGTPWRKRKGIAGPGLPDLSCALVLQPRAQVGTMSPAIALAFLPLVVTLLVRY RHYFRLL VRTVLLRSLRDCLSGLRIEERAFSYVLTHALPGDPGHILTTLDHWSSRCEYLSH MGPVKG QILMRLVEEKAPACVLELGTYCGYSTLLIARALPPGGRLLTVERDPRTAAVAEK LIRLAG FDEHMLVIVGSSSEVIPCLRTQYQLSRADLVLLAHRPRCYLRDLQLLEAHALL PAGATV LADHVLFPGAPRFLQYAKSCGRYRCRLHHTGLPDFPAIKDZIAQLTYAGPG
Source	in vitro E.coli expression system
Target Names	LRTOMT
Protein Names	Recommended name: Transmembrane O-methyltransferase EC= 2.1.1.6 Alternative name(s): Catechol O-methyltransferase 2 Protein LRTOMT2
Expression Region	1-291
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.