



# Recombinant Helicobacter pylori Cbb3-type cytochrome c oxidase subunit CcoP (ccoP)

<b>Product Code</b>	CSB-CF510607HUU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	D0K261
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Helicobacter pylori (strain 52)
<b>Sequence</b>	MDFLNDHINVFLGIAALVILVLTIIYESSSLIKEMRDSKSSQGELMENGHLIDGIGEF ANNV PVGWIASFMCTIVWAFWYFFFGYPLNSFSQIGQYNEEVKAHNQKFEAKWKNL GQKELVDM GGGIFLVHCSQCHGITAEG LHGSAQNLVRWGKEEGIMDTIKHGSKGMDYLAG EMPAMELD EKDAKAIASYVMAEISSVKKTKNPQLIDKGKELFESMGCTGCHGNDGKGLQEN QVFAADL TAYGTENFLRNILTHGKKNIGHMPSFKYKNFSDLQVKALAEFIQSLKPLED
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	ccoP
<b>Protein Names</b>	Recommended name: Cbb3-type cytochrome c oxidase subunit CcoP Short name= Cbb3-Cox subunit CcoP Alternative name(s): C-type cytochrome CcoP Short name= Cyt c(P) Cytochrome c oxidase subunit III
<b>Expression Region</b>	1-292
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.