



# Recombinant Human Methylsterol monooxygenase 1 (MSMO1)

<b>Product Code</b>	CSB-CF620894HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q15800
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Sequence</b>	MATNESVSIFSSASLAVEYVDSLLPENPLQEPFKNAWNYMLNNYTKFQIATWG SLIVHEA LYFLFCLPGFLFQFIPYMKKYKIQKDKPETWENQWKCFKVLNHFHFCIQLPLICG TYYFT EYFNIPYDWERMPRWYFLLARCFGCAVIEDTWHYFLHRLHHKRIYKYIHKVH HEFQAPF GMEAEYAHPLETLILGTGFFIGIVLLCDHVILLWAWVTIRLLETIDVHSGYDIPLN PLNL IPFYAGSRHHDFHMHMNFIGNYASTFTWWDRIFGTDSQYNAYNEKRKKFEKKTE
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	MSMO1
<b>Protein Names</b>	Recommended name: Methylsterol monooxygenase 1 EC= 1.14.13.72 Alternative name(s): C-4 methylsterol oxidase
<b>Expression Region</b>	1-293
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.