



Recombinant Candida albicans Translocation protein SEC62 (SEC62)

Product Code	CSB-CF710595CZD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5AI21
Product Type	Transmembrane Protein
Immunogen Species	Candida albicans (strain SC5314 / ATCC MYA-2876) (Yeast)
Sequence	MSAPTPPQTGQFQVATSAQRSPVAVNIVNYLFQNPILKQRTGLLDNTNDIEFFR FKRLQR ALLSDDYKQKQQNPKNGLIPVANNEDVQKVFVLLIQNQLLLPLQKLHYAEVKAV KGWKPN KEKPTLKRADKAVMDPDVYYGWLYQKPNPYILLYSFLAIAGVFAVILFPLWPNF MKRGVW YLSMGALGLIALFFATAIVRLIYIISLVAFPFPFWLFPNLFEDCGVIESFQPVYAW EEP KKKKGKKKKEPKLEVVEEKVGSSSGDVKVTGSELSNGSSTTGKRKVTLEEVE E
Source	in vitro E.coli expression system
Target Names	SEC62
Protein Names	Recommended name: Translocation protein SEC62
Expression Region	1-293
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.