



Recombinant Mitochondrial prohibitin complex protein 2 (phb-2)

Product Code	CSB-CF344630CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P50093
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MAKQGQEAMKKAIQNARGAGVGLGLVAAAGAAVYGVAQSMFTVEAGHRAIM FNRIGGLST DLYKEGLHFRIPWFQYPIIYDIRARPNQIRSPGSKDLQMVNIGLRVLSRPNPEH LVHIY RTLQGNWEERVLPVICNEVLKGVVAKFNASQLITQRQQVSMMLVRKTLIERALDF NIILDD VSLTELAFSPQYSAAVEAKQVAAQEAQRATFYVERAKQQKQEKIVQAEGEAE SAKLLGEA MKNDPGFLKLRKIRAAQKIARIVSESGNKTYLPTGGLMLNIADTDYLNVTDKRR
Source	in vitro E.coli expression system
Target Names	phb-2
Protein Names	Recommended name: Mitochondrial prohibitin complex protein 2 Short name= Prohibitin-2
Expression Region	1-294
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.