



Recombinant Schizosaccharomyces pombe Protein bem46 (bem46)

Product Code	CSB-CF345731SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P54069
Product Type	Transmembrane Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Sequence	MAGSLSSAIFNVLKYSGMASLAVTLIALGFLYKYQKTLVYPSAFPQGSRENVPT PKEFNM EYERIELRTRDKVTLDSYMLQSESPESRPTLLYFHANAGNMGHRLLPIARVFYS ALNMNV FIISYRGYGKSTGSPSEAGLKIDSQTALEYLMEHPICSKTKIVVYGQSIGGAVAIA LTAK NQDRISALILENTFTSIKDMIPTVFPYGGSIISRFCTEIWSSQDEIRKIKKLPVLFSL GE KDEIVPPPQMVLLFGLCGSAKKKFHSFPKCTHNDTCLGDGYFQVIADFLAENDI NTPAS
Source	in vitro E.coli expression system
Target Names	bem46
Protein Names	Recommended name: Protein bem46
Expression Region	1-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.