



Recombinant *Saccharomyces cerevisiae*

Nucleoporin POM34 (POM34)

Product Code	CSB-CF619602SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q12445
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MKIQAGQLGLDDNDVPGPLPDTDSKPSSQSQNDTPMFKLGNFESPVLKELSR RTV NKEME TQRIMTNVIAFAFWNLLVKFIKFFWNNTHVGRQFCNRLSRIHLYMLTFHTLKKA NIIYHT TFSWLNAELLDYLFHLLISLNILFSLWKLLSTVKVSDLNLTDRQKKLLGVDMQSS VDTGL QPQHPHYVSTSKISQMAQNKTHIPQTNLKNHPAYLFKGLETPLKARQREMAEE QTKLQSQ SLHTKNVFGTLQRHSGISSTLVSANNDNNSPHTPVTRKGYIPSSKYAYMMNSQ SPRGKI
Source	in vitro E.coli expression system
Target Names	POM34
Protein Names	Recommended name: Nucleoporin POM34 Alternative name(s): Nuclear pore protein POM34 Pore membrane protein POM34
Expression Region	1-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.