



Recombinant Drosophila melanogaster Zinc finger protein-like 1 homolog (CG5382)

Product Code	CSB-CF895122DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9VD26
Product Type	Transmembrane Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sequence	MGLCKCPKRLVTNQFCFEHRVNVCEHCMVQSHPKCIVQSYLQWLRDSDYISN CTLCGTTL EQGDCVRLVCYHVFHWDCLNARQAALPANTAPRGHQCPACSVEIFPNANLVS PVADALKS FLSQVNWGRNGLGLALLSEEQSSSLKAIKPKVASQAAVSNMTKVHHIHSGGER ERTKPNG HDAVSPHSVLLMDAFNPPSAGDYASSRRPLLPRQSPIGGTDRDDNKYQRRTP AELFSRWT RRFYAPSSRPPWRRTWFLVTAGILAFVLFVYLMAWLGRGGSDAVDEGWNNP NPQPNHYE
Source	in vitro E.coli expression system
Target Names	CG5382
Protein Names	Recommended name: Zinc finger protein-like 1 homolog
Expression Region	1-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.