



Recombinant UPF0187 protein yneE (yneE)

Product Code	CSB-CF834829EOD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8XAZ3
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O157:H7
Sequence	MIVRPQQHWLRRIFVWHGSVLSKISSRLLLNFLFSIAVIFMLPWYTHLGIKFTLA PFSIL GVAIAIFLGFRRNAGYARYVEARKLWGQLMIASRSLLREVKTTLPDSASVREFA RLQIAF AHCLRMTLRKQPQVEVLAHYLKTEDLQRVLASNSPANRILLIMGEWLAVQRRN GQLSDIL FISLNDRLNDISAVLAGCERIAYTPIPFAYTLILHRTVYLFCIMLPFALVVDLHYMT PFI SVLISYTFISLDCLAELEDPFGTENNDLPLDAICNAIEIDLLQMNDEAEIPAKVLP DRH YQLT
Source	in vitro E.coli expression system
Target Names	yneE
Protein Names	Recommended name: UPF0187 protein yneE
Expression Region	1-304
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.