



Recombinant Walleye dermal sarcoma virus ORF-B protein

Product Code	CSB-CF773330WBH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q88940
Product Type	Transmembrane Protein
Immunogen Species	Walleye dermal sarcoma virus (WDSV)
Sequence	MFSDSDSSDEELSRIITDIDESPQDIQQSLHLRAGVSPEARVPPGGLTQEEWTI DYFSTY HTPLLNP GIPHELTQRCTDYTASILRRASAKMIQWDYVFYLLPRVWIMFPFIARE GLSHL THLLTLTTSVLSATSLVFGWDLTVIELCNEMNIQGVYLPEVIEWLAQFSFLFTHV TLIVV SDGMMDLLLMFPM D IEEQPLAINIALHALQTSYTIMTPILFASPLLRIISCVLYAC GHCP SARMLYAYTIMNRYTGESIAEMHTGFRCFRDQMIAYDMEFTNFLRDLTEEETP VLEITEP EPSPT E
Source	in vitro E.coli expression system
Protein Names	Recommended name: ORF-B protein
Expression Region	1-306
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.