



Recombinant Danio rerio Transmembrane protein 177 (tmem177)

Product Code	CSB-CF023759DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6PBI8
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MSSRFLKISVFIQKYRTPLLLIGCGGVFSANIFYHIFPDHTYKKVYQAWHKGEPA SLSEK LQNIFQEVLKDSSISTSGNFSAFAAFGFHPVGAGVPWLPSGAQIGIPANFNSST ADLEGI TNRTILINGKELEWDSDSGVALKNSLVFSLEAQKFAIAREVARLGSGGPILHAAV APVCL AGACVYSVALKQIFRFQAGSIIFRGVVNLLSLGLGVMTYVLAADSVSQWLDYRS DRRAAG LSRDYAKGGLEFYEKILTRNKTLRSLMGQKGEEMYAPSGNLFPAHLLQLRHAT YTSRRDR ILNLLKNENV
Source	in vitro E.coli expression system
Target Names	tmem177
Protein Names	Recommended name: Transmembrane protein 177
Expression Region	1-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.