



# Recombinant Danio rerio Estradiol 17-beta-dehydrogenase 12-B (hsd17b12b)

<b>Product Code</b>	CSB-CF764830DIL
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q6QA33
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Sequence</b>	MEPFADALFWVGAVTVLWLSVSSLWLSLINGIRVWILGNGNLMRASSLGKWAV VTGATDGI GKAYAEELARRGFAIVLISRTQEKLDEVSKAIESKYKVETKTISADFGSVDIYPKI ESGL AGLEIGVLVNNVGVSYSEFFLNIPDVDSFINNMININIMSVCMTRLVLPRMV DRSKG VILNVASASGMYPVPLLTLYSSTKAFVDFFSRGLDAEYKSKGIIIQSVLPFYVTTK LSKI RKPTLDIPTPERYVKAQLSTIGLQTQSNGYLPHEIMGWVTASLLPAKLLNKYVM GMGLSQ RARYLKKQKQG
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	hsd17b12b
<b>Protein Names</b>	Recommended name: Estradiol 17-beta-dehydrogenase 12-B EC= 1.1.1.62 Alternative name(s): 17-beta-hydroxysteroid dehydrogenase 12-B Short name= 17-beta-HSD 12-B Short name= zf3.3 Short name= zfHSD17B12B
<b>Expression Region</b>	1-311
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.