



Recombinant Mouse NEDD4 family-interacting protein 2 (Ndfip2)

Product Code	CSB-CF856339MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q91ZP6
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	RRSASDAELSAGAEGATGSEAAPPGLGGRTRGGGRGSAAAAATTSTREAE GAERRGDTP ARKPDPEAGRMDHHQLGTGRYQVLHNEEDNSESSAVEQPSTSSLAAPTVEA AASAPALDP DSPPPYSSITVEAPTTSDTDVYSEFYVPPPYSVATSLPTYDEAEKAKAAALAA AAADAP QRNQEEDCTPRDDFSDVEQLRVGNDFMLAFFMAFIFNWLGFCLSFCTNTIA GRYGA CGFGLSLIKWILIVRFSDYFTGYFNGQYWLWWIFLVLGLLLFFRGFVNYLKVRN MSESMA AAHRTRYFFLL
Source	in vitro E.coli expression system
Target Names	Ndfip2
Protein Names	Recommended name: NEDD4 family-interacting protein 2 Alternative name(s): NEDD4 WW domain-binding protein 5A
Expression Region	1-311
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.