



Recombinant Metaxin-1 homolog (mtx-1)

Product Code	CSB-CF528603CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O45503
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MELHIWPSDFGLPTIDVVSLQFLACSKMCASPVRVIQSTRPWRSPSGELPMVA QTEGEAK PVTDFEKFVDILKKCGQDVVIDADLTIEKAQLDAFSCYLHHNLYPAVMHTFWT DELNYN TVTQYWYASHLHFPYNLYYLEKRRRKALRLLAGKNDTEILKEAFMALNTLSTKL GDNKFF CGNKPTSLDALVFGYLAPLLRVPLPNDRLQVQLSACPNLVRVETVSSIYPLG EDELKR QQANRKMWQSRISKAKADKEAAKTTEEASESLPEEPPMRDAILFTLGALTLSL VFAIHTG LIQVSVEEEEISE
Source	in vitro E.coli expression system
Target Names	mtx-1
Protein Names	Recommended name: Metaxin-1 homolog Alternative name(s): Mitochondrial outer membrane import complex protein 1
Expression Region	1-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.