



Recombinant *Thermus thermophilus* Uncharacterized PIN and TRAM-domain containing protein TTHA0540 (TTHA0540)

Product Code	CSB-CF713088TNT
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5SKV3
Product Type	Transmembrane Protein
Immunogen Species	<i>Thermus thermophilus</i> (strain HB8 / ATCC 27634 / DSM 579)
Sequence	DWGLLPQSPSLLSLNRLYLALAGLLTGLLLGPRLEGALEARLKRLRSLPPEVVV ATTLGS TIGLLLAVLLTTLAQQVPGFSPVHSLLLALGLVALFVYLALGYRAYFRLPEPKPAP RGGK VLDTSVLVDGRVAEVAAVGFLEGPLWVPHFVLKELQHFADSQDPLRRAKGRR GLETLELRL REAAPLEVLETPKGESVDEKLLFLARDLEAALVTNDHALLQMARIYGVKALSI QALAQA LRPQLQVGDTLKLLILKEGKEPHQGVGYLEDGSMVVVDGGSRYRGQEIEVVVT QAIQTQV GRLFFARPAQGAQ
Source	in vitro E.coli expression system
Target Names	TTHA0540
Protein Names	Recommended name: Uncharacterized PIN and TRAM-domain containing protein TTHA0540
Expression Region	24-336
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.