



Recombinant Arabidopsis thaliana Uncharacterized mitochondrial cytochrome b-like protein AtMg00590 (AtMg00590)

Product Code	CSB-CF309387DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P93314
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MTIRNQRFSLLKQPISSTLNQHLVDYPTPSNLSYWWGFGPLAGTMILSVLSSPA LVSGLM VARAKNLVHSVLFPPIPIFFSINQLFHYFCRLPIIKHLATKCQLLFLISHLLLLVLT KL VLDLGGYLFMDDLSRALSQFVPGFSGGLGGGSNTPPNPSGDFFLSSYQTSDP DYHDQRRG DSYFSSAPGVQETHRHASGSSTNLHLNLNDQSQDPIFLEVERLSLKCDKVKEK TILKTQS LLERGYHIPDERDIERAINVVMTEHETIDIDRRRKRFFYYLYSCLGKTGNKFWM ELLETL ADYNINIKSDSDN
Source	in vitro E.coli expression system
Target Names	AtMg00590;
Protein Names	Recommended name: Uncharacterized mitochondrial cytochrome b-like protein AtMg00590 Alternative name(s): ORF313
Expression Region	1-313
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.