



Recombinant Human Transmembrane protein 231 (TMEM231)

Product Code	CSB-CF862028HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9H6L2
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MALYELFSHPVERSYRAGLCSKAALFLLAAALTYIPLLVAFRSHGFWLKRSS YEEQPT VRFQHQVLLVALLGPESDGFLAWSTFPAFNRLQGDLRLVPLVSTREEDRNQD GKTDMLHF KLELPLQSTEHVLGVQLILTFSYRLHRMATLVMQSM AFLQSSFPVPGSQLYVN GDLRLQQ KQPLSCGGLDARYNISVINGTSPFAYDYDLTHIVAAYQERNVTTVLNDPNPIWL VGRAAD APFVINAIIRYPVEVISYQPGFWEMVKFAWVQYVSILLIFLWVFERIKIFVFQNQV VTTI PVTVTPRGDLCKEHL
Source	in vitro E.coli expression system
Target Names	TMEM231
Protein Names	Recommended name: Transmembrane protein 231
Expression Region	1-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.