



Recombinant *Saccharomyces cerevisiae* Pheromone-regulated membrane protein 5 (PRM5)

Product Code	CSB-CF461187SVP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B3LTW4
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain RM11-1a) (Baker's yeast)
Sequence	MTVITIAKRGLPKLTTSTSSSTTTASSSSTITSVASSSSSLPLLSNSTSSSIIPSITPP SR NGNPYILDSGDMPNGTVFIVVGGIAGVIFLAILLWWVITTYSSHRLTRSVQDYES KMFST QHTQFYGDSPYMDYPAKENFQDQVHISESDISPGNKDESVKDALVSHTNNEK PFLSNFER PLSSLVSESNRNSLFISPTGDILYKTRLSKLYQESPRLLQKPVIMTSDNVSTNSL VSTIS SSSASSLDNGNEKEVGEDIRKPAKIASSPSRKLNSPESDGSVNRNHSKGNLL VVQSKRK PTPSTYLEHMLEGKEQDE
Source	in vitro E.coli expression system
Target Names	PRM5
Protein Names	Recommended name: Pheromone-regulated membrane protein 5
Expression Region	1-318
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.