



Recombinant *Saccharomyces cerevisiae* Pheromone-regulated membrane protein 5 (PRM5)

Product Code	CSB-CF522245SVK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	E7NIY5
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain FostersO) (Baker's yeast)
Sequence	MTVITIAKRGLPKLTTSTSSSTTTASSSSTITSVASSSSSLPLLSNSTSSSIIPSITPP SR NGNPYILDSGDMPNGTVFIXVGGIAGVIFLAILLWWVITTYSSHRLTRSVQDYES KMFST QHTQFYGDSPYMDYPAKENFQDQVHISESDISPGNKDESVKDALVSHTNNEK PFLSNFER PLSSLVSESNRNSLFISPTGDILYKTRLSKLYQESPRLLQKPVIMTSDNVSTNSL VSTIS SSSASSLDNGNEKEVGEDIRKPAKIASSPSRKLLNSPESDGSVNRNHSGKGNLL VVQSKRK PTPSTYLEHMLEGKEQDE
Source	in vitro E.coli expression system
Target Names	PRM5
Protein Names	Recommended name: Pheromone-regulated membrane protein 5
Expression Region	1-318
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.