



Recombinant GILT-like protein C02D5.2 (C02D5.2)

Product Code	CSB-CF333875CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P34276
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MSKYLRTSHPSPLICRPILTFSSLHILTAFLISGNSSYINWSVTIILLILGLYACHRF LN MKKLTRDAQEPLLPHIRPSTRNKIRKYIYGTFILSIFLLYRSLNSPDTSKHGIGRN GDF IEYEPKAGPTIKEPVENIVKLDVYMEAQCPDTSRFFRQQLKKAWDILGRLNRIEL NVIPF GKARCTEKGNDFECQCQHGPTTECQINQLMNCVIDRFGFPHRYLPGVLCMQG KYSLDEAMK CVTENYPSEYERMRECASGTRGRRLALSGQKTASLTPAIDFIPWIVINGSRNS DALYDL TQNVCEAMQPMP SACKDYLRSLQ
Source	in vitro E.coli expression system
Target Names	C02D5.2
Protein Names	Recommended name: GILT-like protein C02D5.2
Expression Region	1-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.