



Recombinant Friend spleen focus-forming virus Glycoprotein 42 (env)

Product Code	CSB-CF335664FDE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P31793
Product Type	Transmembrane Protein
Immunogen Species	Friend spleen focus-forming virus (strain BB6) (FSFFV)
Sequence	VQRDSPHQVFNVTWRVTNLMTGQTANATSLGTMDFPKLYFDLCDLIGND WDETRLGC RTPGEGKRARTFDLYVCPGHTVPTGCCGPREGYCGKWGCETTQAYWKPS SSWDLISLKR GNTPKDRGPCYDSSVSSGVQGATPGGRCNPLVLKFTDAGKKASWDAPKVV GLRLYRSTGT DPVTRFSLTRQVLNIGPRVPIGPNPVISDQLPPSRPAQIMLPRPPQPPPPGTASI VPETA PPSQQPGRDRLLNLVNKAYQALNLTSPDKTQECWLCLVSRPPYYEGVAVLG TNSNHHTL ISTIMGLLIILLLLLLILLLLWTLHS
Source	in vitro E.coli expression system
Target Names	env
Protein Names	Recommended name: Glycoprotein 42 Short name= gp42
Expression Region	33-356
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.