



Recombinant Uncharacterized protein T25D10.1 (T25D10.1)

Product Code	CSB-CF607330CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q10017
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MFLRRRNLSNSSRIICIISIIVLLLLIIISLYPHKRTQFGRYSTRQKTIRFQHSTGEIGD QL FSLLSHLGVAKTLYRIPVINSANNSKLIDTLNAMFTRFPSILQQFLIAIEPPTAVN REL GIENSSYEDPLTKFSEDTSSSLMVKNGGFKSFKYFDNLRSDIRLWVLEDAESVL EAQNLI TKSQRNNFKICVHATLESNKNCVKAIAQILNHYTNEYEDVMLIIASPFPEFTRFI FTNS RIRKYKTEKFSLISSPEMQIIFSRIYCDVVFLTVPYSTHGWWMGYLAKDDNSH VFYFDP DMFPKNRTANQEDYFPPKWKLSRKIQ
Source	in vitro E.coli expression system
Target Names	T25D10.1
Protein Names	Recommended name: Uncharacterized protein T25D10.1
Expression Region	1-327
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.