



Recombinant Escherichia coli Uncharacterized oxidoreductase YgjR (ygjR)

Product Code	CSB-CF331580ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P42599
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MIRFAVIGTNWITRQFVEAAHESGKYKLTAVYSRSLEQAQHFANDFSVEHLFTS LEAMAE SDAIDAVYIASPNSLHFSQTQLFLSHKINVICKPLASNLAEVDAAIACARENQV VLFEA FKTACLPNFHLLRQALPKVVGKLRKVFFNYCQYSSRYQRYLDGENPNTFNPAPS NGSIMDI GFYCLASAVALFGEPKSVQATASLLASGVDAQGVVMDYGDFSVTLQHSKVS DSVLASEI QGEAGSLVIEKLSECQKVCVPRGSQMQLTQPQHINTMLYEAELFATLVDEH LVDHPGL AVSRITAKLLTEIRRQTGVIFPADSVKL
Source	in vitro E.coli expression system
Target Names	ygjR
Protein Names	Recommended name: Uncharacterized oxidoreductase YgjR EC= 1.-.-.-
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.