



Recombinant Human NEDD4 family-interacting protein 2 (NDFIP2)

Product Code	CSB-CF885735HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9NV92
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MARRRSQRVCASGPSMLNSARGAPELLRGTATNAEVSAAAAGATGSEELPP GDRGCRNGG GRGPAATTSSSTGVAVGAEHGEDSLSRKPDPEPGRMDHHQPGTGTRYQVLLNE EDNSESSAI EQPPTSNPAPQIVQAASSAPALETDSPPPYSSITVEVPTTSDTEVYGEFYPVP PPYSVA TSLPTYDEAEKAKAAAMAAAAAETSQRIQEEECPPRDDFSDADQLRVGNDGIF MLAFFMA FIFNWLGFCLSFICITNTIAGRYGAICGFGLSLIKWILIVRFSDYFTGYFNGQYWL WWIFL VLGLLLFFRGFVNYLKVRNMSESMAAAHRTRYFFLL
Source	in vitro E.coli expression system
Target Names	NDFIP2
Protein Names	Recommended name: NEDD4 family-interacting protein 2 Alternative name(s): NEDD4 WW domain-binding protein 5A Putative MAPK-activating protein PM04/PM05/PM06/PM07 Putative NF-kappa-B-activating protein 413
Expression Region	1-336
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.