



Recombinant DnaJ homolog dnj-2 (dnj-2)

Product Code	CSB-CF615036CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q17433
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MRSAIAAPILFLLVSFFVQECESVGFAPELYCGLENCYDVLEVNREEFDKQKLA KAYRAL ARKHHPDRVKNKEEKLLAEERFRVIATAYETLKDDEAKTNYDYYLDHPDQRFY NYYQYYR LRAAPKVDLRIVIVGTILISLFQFLSAKHKFSEAIEYATGVGKFRNMAIKDGIDKG LLE MDRNGKLLKKNKGVDNDEVIKQIIDNLDVTGGYKRESIYDTLAWHTIIFPLTIFRYI KWT ALWYWRFAIQKEEYDDDAKLYLIRKYIGVSQMEFDQKYTDEDIDDLFERECWL KRNCAW KAERDAAEQEKMAQSGRYKRYKRYMKNAGTISFVDED
Source	in vitro E.coli expression system
Target Names	dnj-2
Protein Names	Recommended name: DnaJ homolog dnj-2 Alternative name(s): DnaJ domain protein 2
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.