



# Recombinant Human herpesvirus 2 Ribonucleoside-diphosphate reductase small chain (UL40)

<b>Product Code</b>	CSB-CF304000HJX
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P69521
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Human herpesvirus 2 (strain HG52) (HHV-2) (Human herpes simplex virus 2)
<b>Sequence</b>	MDPAVSPASTDPLDTHASGAGAAPIPVCPYPERYFYTSQCPDINHLSLSILNR WLETEL VFVGDEEDVSKLSEGELGFYRFLFAFLSAADDLV TENLGGLSGLFEQKDILHYY VEQECI EVVHSRVYNIQLVLFHNNDQARRAYVARTINHPAIRVKVDWLEARVRECDSIP EKFILM ILIEGVFFAASFAAIAYLR TNLLRVTCQSNDLISRDEAVHTTASCYIYNNYLGGH AKPE AARVYRLFREAVDIEIGFIRSQAPT DSSILSPGALAAIENYVRF SADRLGLIHMQ PLYS APAPDASFPLSLMSTDKHTNFFECRSTSYAGAVVNDL
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	RIR2
<b>Protein Names</b>	Recommended name: Ribonucleoside-diphosphate reductase small chain Short name= R2 EC= 1.17.4.1 Alternative name(s): Ribonucleotide reductase 38 kDa subunit Ribonucleotide reductase small subunit
<b>Expression Region</b>	1-337
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.