



Recombinant Spiroplasma virus SpV1-R8A2 B Uncharacterized protein ORF2 (ORF2)

Product Code	CSB-CF325386SJZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P15893
Product Type	Transmembrane Protein
Immunogen Species	Spiroplasma virus SpV1-R8A2 B (SpV1) (Spiroplasma virus 1)
Sequence	MKKFIFFFKNYCYISGSMLLFLSLIDLLLWIISLYCVGLVFWILFVLQCIYFVWWLW KNIF YQLNSFGLVNFVWDNPLSVIIGKLGTKLLTYLSQTMKLLTDEIYSNYPLEDD KVKVL TFKNLDFTDRTPVPPDDSVILFDESYLYIDGTSPHDEKKVHSGKIPWIVLARHF GNRAL FTAQREGMIWNNIRQLASGIIIPISLKKPVIKKGFNFFNRRFFIMQIGIFQDITDYEIW KT KSVERTAEGKRVKHKSDVGLGIRFFKIIIPLEFANKYDSQWLKFVRDLKNDEIVN KKEYY WSEITKLSVKERLELFDIDILKKNLKPKEKNGKDD
Source	in vitro E.coli expression system
Target Names	ORF2
Protein Names	Recommended name: Uncharacterized protein ORF2 Alternative name(s): Gene 2 protein
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.