



# Recombinant *Pyrenophora tritici-repentis* 3-ketoacyl-CoA reductase (PTRG\_11203)

<b>Product Code</b>	CSB-CF451194FHT
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B2WMJ3
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Pyrenophora tritici-repentis</i> (strain Pt-1C-BFP) (Wheat tan spot fungus) ( <i>Drechslera tritici-repentis</i> )
<b>Sequence</b>	MANILEPFGIRIDADNSFVQAAVTGFLLVGIASFAAPLISTIRVLLSLFVLP GKSLT TFG PRGTWALITGASDGIGKEFALS LAAKGYNLILVSRTQSKLDSL SADITSKYGPKI AVKTL AMD FALNKDADYNNM KKLIEGLDVSILINNVGLSHSIPVPFTETPKQEMTDIIMIN CMAT LRVTQLVTPGMVSRKRGLVLTMASFGGFFPTPLLATYSGSKAFLQQWSTALAS ELEPHGV YVQCVQSHLVTTAMSKIRKTSALVPNPKQFVDATLSKIGRSGGAQGV AFTSTP YWSHGLM HWFLSRFLGERSETVVKVNRGMHENIRRRALRKAERDAKKQ
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	PTRG_11203
<b>Protein Names</b>	Recommended name: 3-ketoacyl-CoA reductase Short name= 3-ketoreductase Short name= KAR EC= 1.1.1.- Alternative name(s): Microsomal beta-keto- reductase
<b>Expression Region</b>	1-341
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.