



Recombinant Lipopolysaccharide biosynthesis protein wzzE (wzzE)

Product Code	CSB-CF360350EOD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AG01
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O157:H7
Sequence	MTQPMPGKPAEDAENELDIRGLFRTLWAGKLWIIGMGLAFALIALAYTFFARQE WSSTAI TDRPTVNMLGGYYSQQQLRNLDVRSNMASADQPSVMDEAYKEFVMQLAS WDTRREFWLQ TDYYKQRMVGNKADAALLDEMINSIQFIPGDFTRAVNDSVKLIAETAPDANNL LRQYVA FASQRAASHLNDELKGAWAARTIQMKAQVKRQEEVAKAIYDRRMNSIEQALKI AEQHNIS RSATDVPAAEELPDSEMFLGRPMLQARLENLQAVGPAFDLDYDQNRAMLNTL NVGPTLDP RFQTYRYLRTPEEPVKRDSPPRAFLMIMWGIVGGLIGAGVALTRRCSK
Source	in vitro E.coli expression system
Target Names	wzzE
Protein Names	Recommended name: Lipopolysaccharide biosynthesis protein wzzE
Expression Region	1-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.