



Recombinant *Drosophila erecta* Serine protease HTRA2, mitochondrial (HtrA2)

Product Code	CSB-CF461207DLK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B3P3J9
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila erecta</i> (Fruit fly)
Sequence	AVVSAAVIKREDFTPTIAASKMTGRRRDFNFIADVAGCADSVVYIEIKDTRHFD YFSGQ PITASNGSGFIIQNGLILTNAHVINKPHTMVQVRLSDGRTFPATIEDVDQTSD LATLR IQVNNLSVMRLGKSSTLRSGEWWVALGSPLALSNTVTAGVISSTQRASQELGL RNRDINY LQTDAAITFGNSGGPLVNL DGEAIGVNSMKVTAGISFAIPIDYVKVFLERAAEKR KKGSA YKTGYPVKRYMGITMLTLPDILFELKSRSQNMPSNLTHGVLVWKVIVGSPAHS GGLQPG DIVTHINKKEIKNSSDVYDALADNSKNLDIVILRGVKQMHVITITPEDP
Source	in vitro E.coli expression system
Target Names	HtrA2
Protein Names	Recommended name: Serine protease HTRA2, mitochondrial EC= 3.4.21.108 Alternative name(s): High temperature requirement protein A2
Expression Region	75-422
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.