



# Recombinant *Drosophila pseudoobscura* *pseudoobscura* Serine protease HTRA2, mitochondrial (HtrA2)

<b>Product Code</b>	CSB-CF644931DME
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q297U2
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Drosophila pseudoobscura pseudoobscura</i> (Fruit fly)
<b>Sequence</b>	ALVSAVLIERHRDELTPTISARSLKGRRNEFNFIADVAGCGDSVVYIEIKDTRH FDYFS GQPITASNGSGFVIEQNGLILTNAHVINKPHTMVQVRLSDGRTFPATIEDVDQ TSDLAT LRIQVSGLPVMKLGKSSTLRSGEWWVALGSPLALSNTVTAGVISATQRASQEL GLRNRDI NYLQTDAAITFGNSGGPLVNLGDGEAIGVNSMKVTAGISFAIPIDYVKVFLERAAE RRKKG SAHKTGYPVKRYMGITMLTLTPDILFELKSRSQNMPNNLMHGVLVWKVIVGSP AHSGGLQ PGDIVTHINKKEIKNSSDVYDALAEGRKDLEIVILRGVKQMHSVKITPEDP
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	HtrA2
<b>Protein Names</b>	Recommended name: Serine protease HTRA2, mitochondrial EC= 3.4.21.108 Alternative name(s): High temperature requirement protein A2
<b>Expression Region</b>	78-427
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.