



Recombinant Simian virus 40 Minor capsid protein VP2

Product Code	CSB-CF360955SHP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P03093
Product Type	Transmembrane Protein
Immunogen Species	Simian virus 40 (SV40)
Sequence	GAALTLLGLDIATVSEAAAATGFSVAEIAAGEAAAAIEVQLASVATVEGLTTSEAI AAIG LIPQAYAVISGAPAAIAGFAALLQTVTGVSAVAQVGYRFFSDWDHKVSTVGLYQ QPGMAV DLYRPDDYYDILFPGVQTFVHSVQYLDPRHWGPTLFNAISQAFWRVIQNDIPRL TSQELE RRTQRYLRDSLARFLEETTWTVINAPVNWYNSLQDYYSTLSPIRPTMVRQVAN REGLQIS FGHTYDNIDEADSIQQVTERWEAQSQSPNVQSGEFIEKFEAPGGANQRTAPQ WMLPLLLG LYGSVTSALKAYEDGPNKKRKLRSRGSSQKTKGTSASAKARHKRRNRSSRS
Source	in vitro E.coli expression system
Protein Names	Recommended name: Minor capsid protein VP2 Alternative name(s): Minor structural protein VP2
Expression Region	2-352
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.