



Recombinant Danio rerio Blood vessel epicardial substance (bves)

Product Code	CSB-CF002887DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5PQZ7
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MSNTTSALPSSVPAVSLDPNATLCQDWEQSHHLLFHLANLSLGLGFLIPTTLAL HMIFLR LLLMTGCSLFIAWATLYRCTLDMVWNVVFLLVNFMHFFFLYKRRPIKIDRELK SVYKR MFEPLHVREALFQRLTGQFCTIQTLKKGQVYAAEDKTSVDERLSILLKGKMKVS YRGHFL HNIYTNAFIDSPEFRSTQMNRGERFQVTIAAEENCKLLCWSRERLTYFLESESF LNEVFR YLIGKDITNKLYSLNDPTLSDKAVKKMDRQPSLCSQLSMMQMRNSMASTSDT DDVLNQIL RGGSTGSSLQKNPLTKTSTTMKPIEEGLEDDVFESSEPTTSQNVSKTTKKDI
Source	in vitro E.coli expression system
Target Names	bves
Protein Names	Recommended name: Blood vessel epicardial substance Short name= BVES Alternative name(s): Popeye domain-containing protein 1 Short name= Popeye protein 1
Expression Region	1-352
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.