



Recombinant Escherichia coli O17:K52:H18 UPF0283 membrane protein ycjF (ycjF)

Product Code	CSB-CF487922EOG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7N4C9
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Sequence	MTEPLKPRIDFDGPLEVDQNPKEFRAQQTFDENQAQNFAPATLDEAPEEEGQV EAVMDAAL RPKRSLWRKMVMGGLALFGASVVGQGVQWTMNAWQTQDWVALGGCAAGA LIIGAGVGSVV TEWRRLLWRLRQRAHERDEARDLLHSHGTGKGRAFCCKLAQQAGIDQSHPAL QRWYASIHE TQNDREVVSLEYAHLVQPVLDAAQARREISRSAAESTLMIASPLALVDMAFIAWR NLRLIN RIATLYGIELGYYSRLRFLKLVLLNIAFAGASELVREVGMDWMSQDLAARLSTR AAQGIG AGLLTARLGKAMELCRPLPWIDDDKPRLGDFRRQLIGQVKETLQKGKTPSEK
Source	in vitro E.coli expression system
Target Names	ycjF
Protein Names	Recommended name: UPF0283 membrane protein ycjF
Expression Region	1-353
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.