



Recombinant Human GRAM domain-containing protein 2 (GRAMD2)

Product Code	CSB-CF818231HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8IU Y3
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MTALSRSEATEEGGNQQMHRKTASLNSPVSCKEKPDRVEEPPDYSLHWPEG LKGEEIKKC GREGITLNKYNQQYHKLFKDVPLEEVVLKVCSCALQRDFLLQGRLYISPNWLC FHASLFG KDIKVVIPVSVQMIKKHKMARLLPNGLAITTNTSQKYIFVSLLSRDSVYDLLRR VCTHL QPSSKKSLSVREFSGEPESLEVLIPEMKWRKVCPSRSLSLPDNIPCI PPSSVD STDSFF PSRKPPMSEKSRAQVASENGGRWAWPMPGWGPACPKMPNCSPTAKNAV YEEDELEEEPR STGELRLWDYRLLKVFFVLICFLVMSSSYLAFRISRLEQQQLCSLSWDDPVPGH R
Source	in vitro E.coli expression system
Target Names	GRAMD2A
Protein Names	Recommended name: GRAM domain-containing protein 2
Expression Region	1-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.