



Recombinant Pig Protein ATP1B4 (ATP1B4)

Product Code	CSB-CF859782PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9BDK6
Product Type	Transmembrane Protein
Immunogen Species	Sus scrofa (Pig)
Sequence	MRRQLRSRRAPAFPYGYGYRLDDQDEVNQNylaDEEEEEAEEARVMVVPDLE EEEEEEEE KEEDEKEEEESHHQDTRSAAWWQKLQIVNEYLWDPEKRMSLARTGQSLSLLLV IYFFFYAS LAAVITLCMYTLFLTISPYPVPTFTERVKPPGVMIRPFAHSLNFNFNVSEPTWQ HYVISL NGFLQGYNDLQEEEMNVDCPPGQYFIQDGEDDEDKKACQFKRSFLKNC SGL EDPTFGYST GQPCILLKMNRIVGFRPELGDPVKV SCKVQRGDENDIR SISYYPESASFDLRY Y PYYGKL THVNYTSPLVAMHFTDVVKNQAVPVQCQLKGGKIINDVINDRFVGRVIFTLNIET
Source	in vitro E.coli expression system
Target Names	ATP1B4
Protein Names	Recommended name: Protein ATP1B4 Alternative name(s): X,K-ATPase subunit beta-m X/potassium-transporting ATPase subunit beta-m
Expression Region	1-355
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.