



Recombinant Arabidopsis thaliana RING-H2 finger protein ATL51 (ATL51)

Product Code	CSB-CF882935DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9SRQ8
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MGSTGNPNPWGTTYDSYRDCSQGVCSVYCPQWCYVIFPPPPSFYLDDEDSS SSSDFSPLL IALIGILASAFILVSYITLISKYCHRRRHNSSTSAAINRISSDYTWQGTNNNNN NGAT NPNQTIGGGGGDGLDESLIKSITVYKYRKMDGFVESSDCSVCLSEFQENESLR LLPKCNH AFHVPCIDTWLKSHSNCPLCRAFIVTSSAVEIVDLTNQQIVTENNSISTGDDSVV VVNLD LENSRSRNETVNEGSTPKPEMQDSRDGEERRSASLNSGGGVVSIADILREIE DDEESAG VGTSRWVEEGEGEKTPPPSGSAANQTNGISNFLVRSSMAAMKRSGYDRAKN YRLPK
Source	in vitro E.coli expression system
Target Names	ATL51
Protein Names	Recommended name: RING-H2 finger protein ATL51
Expression Region	1-356
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.