



Recombinant Arabidopsis thaliana UPF0496 protein At3g57100 (At3g57100)

Product Code	CSB-CF862964DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9M1J0
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MKFCNPFLRLILRRKPISMISPTMFLSSPTTSSSIGHSNGMLRSRSEDNLNKLG ECMDN MDEDVSKLFIQQTDLREDPDFRLLNHFFTSKGVSQLCESLRTCLERSEN NECLLD EALVDFELEKLGYGGSLEEASFRKTYRDLRNFNAFYNNNSGEEDLDYCEFLR KFQTCHE ELAKMVVKLEKTMKDIDKKLRRVRGRRAIVTAALLAPVIAVIFLSKLVAGLVPIEG LSTF VASRWKSTESLKREKTAMSSMERGIIVALKQVEKISKLVSRLESVERSISLTAE FAVKK RSSVVVAMREVEEERKRLKSTLVLDRETGLCNGFAQFGRTVALEKITEFLSR GDKSSK
Source	in vitro E.coli expression system
Target Names	At3g57100
Protein Names	Recommended name: UPF0496 protein At3g57100
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.