



Recombinant Human Peroxisome assembly protein 12 (PEX12)

Product Code	CSB-CF017798HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O00623
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAEHGAHFTAASVADDQPSIFEVVAQDSLMTAVRPALQHVVKVLAESENPTHY GFLWRWFD EIFTLLDLLLQQHYLSRTSASFSENFYGLKRIVMGDTHKSQRLASAGLPKQQLW KSIMFL VLLPYLKVKLEKLVSSLREEDEYSIHPPSSRWKRFYRAFLAAYPFVNMAWEGW FLVQQLR YILGKAQHHSPLLRLAGVQLGRITVQDIQALEHKPAKASMMQQPARSVSEKIN SALKKAV GGVALSLSTGLSVGVFFLQFLDWWYSSSENQETIKSLTALPTPPPPVHLDYNSD SPLLPKM KTVCPICRKRTRVNDTVLATSGYVFCYRCVFHYVRSHQACPITGYPTVEVQHILKL YSPEN
Source	in vitro E.coli expression system
Target Names	PEX12
Protein Names	Recommended name: Peroxisome assembly protein 12 Alternative name(s): Peroxin-12 Peroxisome assembly factor 3 Short name= PAF-3
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.