



Recombinant Pongo abelii Cell cycle control protein 50A (TMEM30A)

Product Code	CSB-CF732633PYX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5R6C0
Product Type	Transmembrane Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Sequence	AMNYNAKDEVDGGPPCAPGGSAKTRRPDNTAFKQQRLPAWQPILTAGTVLPI FFIIGLIF IPIGIGIFVTSNNIREIEIDYTGTEPSSPCNKCLSPDVTPCICTINFLEKSFEGNVF MY YGLSNFYQNHRRYVKSRDDSQLNGDSSALLNPSKECEPYRRNEDKPIAPCGAI ANSMFND TLELFLIGNDSYPIPIALKKKGIWWTDKNVKFRNPPGGDNLKERFKGTTKPVN WLKPVY MLDSDPDNNGFINEDFIVWMRTAALPTFRKLYRLIERKSDLHPTLPAGRYSLNV TYNYPV HYFDGRKRMILSTISWMGGKNPFLGIAYIAVGSISFLLGVVLLVINHKYRNSNT ADITI
Source	in vitro E.coli expression system
Target Names	TMEM30A
Protein Names	Recommended name: Cell cycle control protein 50A Alternative name(s): Transmembrane protein 30A
Expression Region	2-361
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.