



Recombinant Enterobacteria phage I2-2 Gene 1 protein (I)

Product Code	CSB-CF324467ECP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P15418
Product Type	Transmembrane Protein
Immunogen Species	Enterobacteria phage I2-2 (Bacteriophage I2-2)
Sequence	<p>MAVYVVVTGKLGAGKTLVAVSRIQRTLAKGGIVATNLNLKLHHPQVGRYAKQC RVMRIAD KPTLEDLEAIGRGNLSYDESKNGLIVLDECGTWFNSRNWSDKSRQPVIDWFLH ARKLGWD VIFIIQDISLMDKQAREALAEHVVYCRRLDKLNIPPIIGGLISVLSGGRLPLPKVHFG IVK YGDNPQSLTVDKWIYTGTDLYAAYDTKQIFTSRELSPFPCVSPYYTHGIFAV KRDAKY YMRMTKIYFKKMNRVWLMASFLALGAGVGGFFYKSRQINEQLSNMPVASAQAN TTKTDHTI DELPRLSINSFAQMGYDVNVSFKDAKGKIYYSFDLMKSGYALDIKDSCHITLRK RNYIQQ VTCEG</p>
Source	in vitro E.coli expression system
Target Names	I
Protein Names	Recommended name: Gene 1 protein Alternative name(s): G1P
Expression Region	1-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.