



Recombinant Human E3 ubiquitin-protein ligase RNF149 (RNF149)

Product Code	CSB-CF822755HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NC42
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	RALEWFSAVVNIEYVDPQTNLTVWSVSESGRFGDSSPKEGAHGLVGVPWAP GGDLEGCAP DTRFFVPEPGGRGAAPWVALVARGGCTFKDKVLVAARRNASAVVLYNEERY GNITLPMSH AGTGNIVVIMISYPKGREILELVQKGIPVTMTIGVGTRHVQEFISGQSVVFVAIAFI TMM IISLAWLIFYIQRFLYTGSQIGSQSHRKETKKVIGQLLLHTVKHGEKGIDVDAEN CAVC IENFKVKDIIRILPCKHIFHRICDPWLLDHRTCPMCKLDVIKALGYWGEPGDVQE MPAP ESPPGRDPAANLSLALPDDDGSDDSSPPSASPAESEPQCDPSFKGDAGENTA LLEAGRSD SRHGGPIS
Source	in vitro E.coli expression system
Target Names	RNF149
Protein Names	Recommended name: E3 ubiquitin-protein ligase RNF149 EC= 6.3.2.- Alternative name(s): DNA polymerase-transactivated protein 2 RING finger protein 149
Expression Region	33-400
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.