



Recombinant Mouse Protein FAM70B (Fam70b)

Product Code	CSB-CF682909MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5FW56
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MVGAEGRWREKGNDRGRKRRGGVEESGHRAAFLETAGTLKPNHRAGGA SGHCPPSSGGH FGNAAAPTARARPSGSAGHHRGLCQEEEDLSVVCGLSPAGVYSHPDHWPC NHQDREDCD GRLLPRDYFGAAIVFLSFGVVAAFCCAIVDSVFAARHIEPRPLSAGRCQFYSSG AGYLYD VYQTEVTCYSLNGRCQLKVRNNTCYCCDLYACGSTEPSPAYYEFVGVRSQD VVHLYRLL WVSTVLNVLGLCLGIVTAAVLGAFKDMVPLSQLAYGPSPPPQILYNPAQQILAY TGLCPP SMGVPTCSSYPLPLQPSSAPPASASADLSLPEDTESPSQCQPSRGCSHAPSP CTPAYFLP GEKPPPYAP
Source	in vitro E.coli expression system
Target Names	Tmem255b
Protein Names	Recommended name: Protein FAM70B
Expression Region	1-369
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.